



新一代设备

# classic-cart

新产品

2024



## 技术细节

### 负载能力

总负载能力	150 kg / 330 lbs
层架的承重能力	50 kg / 110 lbs 分布式负载

### 尺寸 - 基本框架

宽度(毫米)	663 / 763 / 1263
深度(毫米)	660
总高度(毫米)	840 / 968 / 1128 / 1288 / 1448 / 1608 / 1768

### 尺寸 - 层架

宽度(毫米)	500 / 600 / 1100
深度(毫米)	491

### 电气化

隔离变压器	1200 VA / 1600 VA / 2000 VA
总开关	yes
插排	yes

Get in touch >>

Europe: +49 89 614425-0 | sales@itd-cart.com  
 America: +1 800 947-3901 | salesusa@itd-cart.com  
 China: +86 21 5031 8818 | saleschina@itd-cart.com  
 India: +91 86 0015 4064 | salesindia@itd-cart.com  
 Australia: +61 3 9330 5050 | salesaustralia@itd-cart.com

www.itd-cart.com



DS.0270.799 - 10/2024. 德国印刷。保留所有修改权。

## 定制设计

### 新颜色

#### 基本色



RAL 9016  
traffic white

#### 装饰条颜色



RAL Design 2507025  
meadow blue



RAL 5017  
traffic blue



RAL 7016  
anthracite grey

### 新设计

个性化企业设计



# 功能优化的高品质设计



**新护角 -**  
设计优雅, 卫生性能更佳



**推动手柄 -**  
手感舒适, 可定制符合  
企业品牌标准

**平滑的单立柱 -**  
具有出色的卫生功能  
和独特的设计



**定制设计 -**  
识别度高

**新的支臂位置 -**  
活动范围大



**额外的背部把手 -**  
便于从任何位置  
对推车进行调整

**后门 -**  
便于存取设备和电缆,  
隐蔽式电缆管理

**国际插座条 -**  
畅通无阻



**配重 -**  
隐蔽式安装, 灵活与隔离变压器相结合

