

pro-cart

Planning Support

Equipment

Castor diameter: D _____ mm

Conductivity: yes no

Isolating transformers: 660 VA with 4-way EURO socket strip
 1200 VA with 6-way EURO socket strip

Insulating monitor: yes no


Socket strip: ZV. _____

Accessories: _____

Design

Base colour: RAL 7035 (Light grey)

Trim colour: .903 (RAL 7035 Light grey)
 .908 (RAL 5015 Sky blue)
 .909 (RAL 1003 Signal yellow)



Dimensions

Total height TH:

Units (U), 1 U = 32 mm	Castors D 100 mm	Castors D 125 mm
21 U	884 mm	909 mm
30 U	1172 mm	1197 mm
35 U	1332 mm	1357 mm

Total width: 610 mm

Total depth: 620 mm

Project _____

Company _____

Street _____

City _____

Fon _____

Fax _____

E-Mail _____

Processor _____

Please note:

Just draw in your selection of accessories at your choice of position between the two dotted lines. Accessories like shelves, drawers etc. for **pro-cart** can be assembled starting **at the first unit (at 32 mm)** from the base.

iTD GmbH
 Lilienthalstr. 19
 85579 Neubiberg, Germany
 Tel. +49 89 61 44 25 100
 E-Mail: sales@itd-cart.com
 Web: www.itd-cart.com



pro-cart

electric height adjuster
 mechanical height adjuster

Equipment


Supply / loading unit: Input voltage:
 AC 110 V / 230 V, 50 / 60 Hz
 Output voltage: DC12 V – 24 V

Accessories: Shelf (PT. _____)
 Block of drawers (PT. _____)

Design

Base colour: RAL 7035 (Light grey)

Trim colour: .903 (RAL 7035 Light grey)
 .908 (RAL 5015 Sky blue)
 .909 (RAL 1003 Signal yellow)



Dimensions

Total height TH:	1357 mm
Height adjuster:	1276 mm to 1576 mm
Installation height:	35 U (1 U = 32 mm)
Total width:	610 mm
Total depth:	620 mm
Castors:	4 twin swivel castors, D 125 mm, 2 with wheel locks, and wall protectors

Project _____

Company _____

Street _____

City _____

Fon _____

Fax _____

E-Mail _____

Processor _____

iTD GmbH
 Lilienthalstr. 19
 85579 Neubiberg Germany Tel.
 +49 89 61 44 25 100
 E-Mail: sales@itd-cart.com
 Web: www.itd-cart.com



Please note:

Just draw in your selection of accessories at your choice of position between the two dotted lines. Accessories like shelves, drawers etc. for **pro-cart with power** supply can be assembled starting **at the fifth unit (at 160 mm)** from the base.